

# Acleara™

## Acne Clearing System



### Visible Improvement:

- Reduction in peri-lesional erythema and a visible calming of lesion

### FDA Cleared for Most Types of Mild to Moderate Acne:

- Mild to Moderate Inflammatory Acne (Acne Vulgaris)
- Comedonal Acne
- Pustular Acne
- Cleared for use on all skin types

### Hygienic and Comfortable Treatments:

- Reduced risk of cross contamination
- No anesthetics, numbing or cooling gels required

### Fast Treatments:

- Full face treatment in under 15 minutes

### Highly Economical:

- Eliminates external cool sprays and anesthetics costs
- Low cost of ownership



BEFORE



AFTER

1 WEEK POST ONE TREATMENT

- Advanced Filtered Broadband Light



# Acleara Mechanism of Action

The Acleara System is powered by advanced vacuum and filtered broadband light technology to treat mild to moderate acne (including acne vulgaris), comedonal and pustular acne.



BEFORE



AFTER



BEFORE



AFTER

BOTH PHOTOS:  
1 WEEK POST ONE TREATMENT

## Position



Acleara handpiece makes contact with treatment area

## Cleanse



Vacuum deep cleans pores by extracting build up of sebaceous material

## Treat



Along with targeted heating of sebaceous glands, the endogenous light activates porphyrins to destroy P. acnes and reduces sebum production through photodynamic action.

## Complete



Obstruction in the pilosebaceous structures are removed and the skin returns to normal.

# ACLEARA™ TECHNICAL SPECIFICATIONS

Fluence	Up to 25 J/cm <sup>2</sup>
Spot Size	13 mm x 25 mm

©2010 Palomar Medical Technologies, Inc. Palomar and Acleara are registered trademarks of Palomar Medical Technologies, Inc. ALL RIGHTS RESERVED.